

Diabetes & **HbA1c** testing



 **- A1c H5**

HbA1c Analyzer (HPLC)




HbA1c Analyzer (HPLC)



Q-A1cH5 fully automated HbA1c Analyzer supports fast HbA1c results output in 130 seconds without Hb variant interference. It provides the outstanding solution for reliable diabetic monitoring. No sample preparation and walk-away operation.

The Elements of Hemoglobin



Gold Standard of Diabetes Diagnose

Glycosylated hemoglobin (HbA1c) is widely recognized as a Gold Standard to monitor diabetes, which can indicate the average plasma glucose concentration over 8 ~ 12 weeks.

HPLC Methodology

High-Pressure Liquid Chromatography (HPLC), to separate HbA1c directly with measuring the absorbance points continually to form chromatogram. Using normal distribution curve fitting auto-iterative algorithm to get precise HbA1c testing result, excluding interference of variant and unstable hemoglobin like LA1c. Standard Analysis Mode will report HbA1a, HbA1b, HbF, LA1c, HbA1c, HbA0 peak areas and percentage. And the result also includes IFCC, NGSP and ADAG value for diverse client needs.



HPLC Technology – Gold Standard Methodology

- HbA1c Results within 130 Seconds
- 5 samples loading capacity, suitable for medium to small labs

Fully Automated - To Minimize Operation Steps


- No sample preparation
- Fully automated system cleaning after test
- External barcode scanner for sample identification

Precise and Reliable – To Serve You Consistently

- Inter measuring $CV \leq 1.5\%$ & Intra measuring $CV's \leq 3\%$
- Superior quality chromatographic resolution to eliminate interferences

Dual Wavelength Detection – To Avoid Interference

- To avoid the reagent peak interference
- More anti-interference abilities, the mutation factor interference to the peak can be easily counteracted
- To eliminate the nonspecific absorption of hemoglobin




Degasser – For Better Result Accuracy

- More Stable Pressure, More Accurate Flow Rate
- To Reduce Background Absorption and Improve Detection Sensitivity
- To Improve the Separation Effect of Column and Prolong Its Lifetime

Compact Size – To Minimize Space Requirements

- Small Footprint Reduces Bench Space Needed
- Most compact fully automated HPLC system



Correlation between A1cH5 and D-10



HbA1c Analyzer (HPLC)

Technical Specifications

Methodology	High - Performance Liquid Chromatography (HPLC)
Test Range	3% - 18%
Precision	CV \leq 1.5%
Test Speed	130 Secs / Test
Sample Type	Venous Blood, Finger Peripheral Blood, Lyophilized Whole Blood 10 μ L(whole blood), 400 μ L (Diluted blood)
Auto Sample Station	5 Positions
Photometer	415nm + 500nm Detector
Filter	\geq 800 Tests
Display	10.1 " TFT True Color LCD Touch Screen
Software	Linux Software with Self - Diagnosis to Monitor and Detect System Errors
QC Curve	Yes, Available for two level control
Reagent Kit	Eluent A, Eluent B, Hemolysin L, Calibrator; QC Materials
Information Input	Scanner or Touch Keypad
Storage	4000 Sample Results
Connection	USB, LAN, LIS Compatible
Printer	Inbuilt Thermal Printer & External Laser Printer
Operation	Temperature 10 ~ 30 °C (50 ~ 86 °F)
Humidity	\leq 80%
Power	AC 100-240V 50/60HZ 120VA
Dimensions	450mm(L)x360mm(W)x540mm(H) (17.7inchx14.17inchx21.2inch)
Weight	32.8KG

QH5-V1.0.06012025

This brochure may be revised or replaced by Q-Line Biotech at any time without notice

Q-LINE BIOTECH PRIVATE LIMITED

C -108, Naraina Industrial Area, Phase-1, New Delhi - 110028, INDIA.

Phone No.: 011 4557 7407 | Toll Free No.: 1800 123 0079

Email: sales@qlinebiotech.com | Website: www.qlinebiotech.com



FOLLOW US

