

Celltac α +

Automated Hematology and ESR Analyzer

MEK-1305

NIHON KOHDEN

Celltac α +



**3-part differential
hematology analyzer
with ESR**



Fighting Disease with Electronics

NIHON KOHDEN

Integration


Transforming the possibility of IVD solutions

Infectious disease is spreading all over the world and becoming a serious problem. For example, 10 million people developed Tuberculosis (TB) which is one of the most serious infectious diseases, and TB caused an estimated 1.6 million deaths in 2017 according to data from the WHO.

What are your laboratory challenges?

What is required for a better clinical outcome?


Celltac α +, equipped with ESR, can help to achieve a better clinical outcome.



CiRHEX (Cell counter integrated rheometric excellence) technology can provide ESR result highly correlated with Westergren method by using HCT value and MCV value from CBC measurement and also RBC aggregation phenomenon.



Nihon Kohden's unique CiRHEX Technology helps you to achieve a better clinical outcome



Result displayed in 2 minutes by single EDTA tube for both CBC and ESR

A single EDTA tube can be used for both CBC and ESR measurement on our Celltac α+, and you will be able to get CBC results on the screen in 1 minute, and an ESR result in 2 minutes, with a single aspiration. This leads to reducing your workload, avoiding the risk of infection and providing a quick report to the patient.



Only 80 µL sample volume required for both CBC and ESR

Unlike the traditional methods for ESR testing, our Celltac α+ requires only 80 µL of blood sample for both CBC and ESR measurement. Not only does this small blood collection volume improve your laboratory workflow, it also improves the patient experience.



No additional reagents, no additional cost

As it is based on conventional hematology analyzers, diluent, detergent and hemolysing reagent are used for CBC and 3 part differential measurement on Celltac α+, but no additional reagents are required for ESR measurement.




Innovation

Maximizes laboratory productivity

Quality hematology testing

DynaHelix Flow technology perfectly aligns RBC and PLT cells for high impedance counting precision using an advanced hydrodynamic-focused sheath flow before passing through the aperture. In addition, the DynaHelix Flow totally prevents the risk of coincidence or re-entry of counted blood cells into the aperture, using the unique DynaHelix Flow stream.

This newly-developed advanced DynaHelix Flow Technology greatly improves counting precision and accuracy.



Integrated QC program



- QC program for ESR is available
- The same QC material can be used for CBC, 3 part diff and ESR
- QC lot management up to 25
- Assay value registration using a handy barcode reader (standard accessory)
- Automated judgement function (pass or fail)
- QC management by assay value, average value or Westgard multirule
- QC graph display and printout (optional)
- Automated calculation of statistical information such as average and SD

Reagent management



Standard accessory, barcode reader



Celltac α+ reagent management system helps easier reagent bottle management with a unique barcode labeled on each reagent. Through this system and use of genuine Nihon Kohden reagents, testing quality is always maintained at a high level.

Thirty-one measuring parameters including ESR and other research parameters

Traditional CBC parameters, WBC 3 part differential parameters, Mentzer Index and RDWI, which are considered to be useful for Thalassemia screening. ESR and other parameters which are related to ESR are available on Celltac α+.

Smart
ColorRun Assist

YELLOW

A panic value (far outside the normal range) needs to be reported to a doctor immediately

ORANGE

Possibly incorrect data due to problems caused by the state of the blood sample or the measuring procedure

RED

Possibly incorrect data due to a technical problem with the instrument or measuring procedure

The figure displays three screenshots of a medical device interface, each showing a patient's lab results. The interface is divided into sections for patient information, sample information, and a list of test results. The results are color-coded to indicate the status of the measurement:

- YELLOW:** A panic value (far outside the normal range) needs to be reported to a doctor immediately. In the screenshot, the MCHC result is 22.0 g/dL, which is highlighted in yellow.
- ORANGE:** Possibly incorrect data due to problems caused by the state of the blood sample or the measuring procedure. In the screenshot, the MCHC result is 22.0 g/dL, which is highlighted in orange.
- RED:** Possibly incorrect data due to a technical problem with the instrument or measuring procedure. In the screenshot, the MCHC result is 22.0 g/dL, which is highlighted in red.

* ASTM protocol is available

```
graph LR; RBC[Red Blood Cell Flags] --- RBC_List["■ Erythrocytosis  
■ Anaemia  
■ Anisocytosis  
■ Microcytosis  
■ Macrocytosis  
■ Hypochromia  
■ Abnormal MCHC"]; WBC[White Blood Cell Flags] --- WBC_List["■ Leucocytosis  
■ Leukopenia  
■ RBC Ghost  
■ Cannot Analyze LY  
■ Cannot Analyze GR  
■ Cannot Differentiate WBC  
■ Large Nucleated Cells  
■ Poor Hemolization"]; PLT[Platelets Flags] --- PLT_List["■ Thrombocytosis  
■ Thrombocytopenia  
■ Platelet Clumps  
■ PLT-RBC Interference"]; ESR[ESR Flags] --- ESR_List["■ High ESR  
■ Increased RBC Aggregation  
■ Decreased Aggregation"];
```

Red Blood Cell Flags

- Erythrocytosis
- Anaemia
- Anisocytosis
- Microcytosis
- Macrocytosis
- Hypochromia
- Abnormal MCHC

White Blood Cell Flags

- Leucocytosis
- Leukopenia
- RBC Ghost
- Cannot Analyze LY
- Cannot Analyze GR
- Cannot Differentiate WBC
- Large Nucleated Cells
- Poor Hemolization

Platelets Flags

- Thrombocytosis
- Thrombocytopenia
- Platelet Clumps
- PLT-RBC Interference

ESR Flags

- High ESR
- Increased RBC Aggregation
- Decreased Aggregation

