


Complete Advantage of Column Agglutination Technology


- Reads*** ▶
- Interprets*** ▶
- Records*** ▶




Documentation of end results is one of the advantage for which Gel Technology is widely accepted and used. Tulip Introduces Matrix™ Card Reader for automatic Reading, Interpretation and Documentation of results. Matrix™ Card Reader is operated through unique Matrix™ Pro-Scan software. Matrix™ Gel Cards are read and identified through unique Barcode on the cards. Matrix™ Card Reader and Matrix™ Pro-Scan along with wide range of Matrix™ Gel Cards are perfect combination for low to medium work load Blood Banks. Combination of Matrix™ Octopus Gel Cards and Matrix™ Card Reader makes it ideal for Diagnostic Laboratories.

Salient Features of Matrix™ Card Reader and Matrix™ Pro-Scan


- ▶ Simple and user friendly software.
- ▶ Compatible with 6 and 8 microtube Matrix™ Gel Cards.
- ▶ Capable of reading Strength of Reaction and Interpretation of Results.
- ▶ Positive identification of Matrix™ Gel Cards through Barcodes.
- ▶ Complete traceability of Test Results, Sample and Operator.
- ▶ Storage and printing of test results along with original image.
- ▶ LIS Compatible.




▶ Card barcode is read by the system.




▶ Easy test selection especially for the cards where different tests can done in one card.



▶ Simple Patient data entry. One window for all samples.



▶ Results window for Blood group test, with grading of reaction and final result interpretation.



▶ Results window for Matrix™ Octopus Complete Grouping Card with grading of reaction and final result interpretation.



▶ Results print preview for crossmatch test, with grading of reaction and final result interpretation.

*Specifications are subject to change without prior notice.

▶ Small ▶ Simple ▶ Efficient

ISO 13485
Certified Company



TULIP DIAGNOSTICS (P) LTD.

For further information contact :

Gitanjali, Tulip Block, Dr. Antonio do Rego Bagh, Alto Santacruz, Bambolim Complex Post Office, Goa - 403 202, INDIA.
Tel.: +91 832 2458546-50 Fax : +91 832 2458544 E-mail : sales@tulipgroup.com Website : www.tulipgroup.com