



matrixTM
GEL•SYSTEM

Card Warmer

CW-2400

KEY FEATURES

- Seven segment LED display with front panel keys.
- Working temperature 37°C (set temp.)
- Warm up time maximum 30 minutes.
- Timer with alarm (15 minutes).

TECHNICAL SPECIFICATIONS

| | |
|-----------------------------|---------------------------------|
| Principle | Incubation |
| Card capacity | 24 cards |
| Operating conditions | Temperature- Ambient 25° - 40°C |
| Power Supply | 230V, 50 Hz. |
| Dimensions | 29 cms x 25.5 cms x 14.5 cms |
| Net Weight | 5.1 Kgs |
| Cat. Nos. | 821CW002400 |

LIST OF STANDARD COMPONENTS SUPPLIED ALONG WITH THE INSTRUMENT

| Components | Qty |
|-------------------|------------|
| Power Cord | 1 No. |
| User Manual | 1 No. |