


# A1Chek

Smart HbA1C Meter



## Certification



\*\*Manufactured in USA FDA approved facility



Easy HbA1c Testing



Accurate results within 5 minutes



Automatic Calibration



No Refrigeration required



Wide operating temperature range of 10°-40°C



Voice Assisted Guidance



Bluetooth Connectivity

**Biosense**  
A Tulip Diagnostics Company

# SIMPLE CONNECTED PRECISE

# A1Chek

Smart HbA1C Meter



Features	Specifications
Test Principle	Boronate Affinity Chromatography
Parameter	HbA1c, eAG (estimated average glucose)
Measuring Range	4 to 14%
Sample Type	Capillary or Venous Blood
Sample Volume	3 $\mu$ L
Precision	CV <3%
Memory	1000 tests
Test Time	5 minutes
Operating Temperature	10° to 40° Centigrade
Certification	NGSP, IFCC, ISO, CE
Power Source	3.7 V lithium ion rechargeable battery
Weight	99.7 gm without battery
Connectivity	Bluetooth
Audio Assistance	Yes

\*Specifications are subject to change without prior information.

Manufactured by :

**Biosense Technologies Pvt. Ltd.**

Bldg no A9, Unit No. 39,40, Pimpalas Village, Mumbai Nashik Highway, Pimpalas, Thane, Maharashtra (India) - 421302

Marketed by :

**TULIP DIAGNOSTICS (P) LTD.**

Gitanjali, Tulip Block, Dr. Antonio do Rego Bagh, Alto Santacruz, Bambolim Complex Post Office, Goa - 403 202, INDIA.

Tel.: +91 832 2458546-50 Fax : +91 832 2458544 E-mail : sales@tulipgroup.com Website : www.tulipgroup.com



**Biosense**  
A Tulip Diagnostics Company