



ERYBANK[®]	ANTI-A₁ LECTIN
Application	For differentiating A ₁ and A ₂ subgroups
Titre	<ul style="list-style-type: none"> • ~ 1:32 with A₁ cells • ~ 1:16 with A₁B cells
Class of Agglutinin	Phytohaemagglutinin
Source	Dolichos biflorus seed extract
Specificity	A ₁ antigen
Standardisation	Follows AABB and FDA guidelines
Storage / Stability	2-8°C, 24 months

Presentation	Pack	Cat.No.
ANTI-A₁ Lectin reagent	5 ml	10200005
	10 ml	10200010

