



ERYCARD™ 2.0 Blood Grouping Card Test for Forward Grouping with Autocontrol	
Application	<ul style="list-style-type: none"> • Bed side ABO & Rh testing of Patient and Donor's blood unit before transfusion • Donor's blood group screening in outdoor donation camps • Outdoor blood grouping camps • Routine laboratory testing • Preconfirmation test, especially at nursing homes • For releasing compatible blood for emergency transfusions in <ul style="list-style-type: none"> • Defence Camps • Ambulance • Casualty services
Principle	Self performing, based on the principle of lateral flow guided by capillary action.
Sensitivity & Specificity	100% correlation with classical tube technique
Reagents	Reagents are pre-dried at the appropriate sample pad Anti-A (IgM) monoclonal antibodies in sample well A Anti-B (IgM) monoclonal antibodies in sample well B Anti-D (IgM) monoclonal antibodies in sample well D Autocontrol without any antibodies (Ctrl) Reagent Buffer (Proprietary Buffer)
Standardisation	Follows FDA guidelines
Test duration	Less than five minutes
Storage / Stability	2-30°C, 24 months

Presentation	Pack	Cat.No.
ERYCARD™ 2.0 kit contains	24 Cards	102800024

Individually pouched device		
Reagent buffer		

