



	Rapid Test for detection of antibodies to <i>Helicobacter pylori</i> in human serum/plasma
Application	For the detection of antibodies to <i>Helicobacter pylori</i> in human serum or plasma.
Principle	Self performing, two site sandwich immunoassay, lateral flow device.
Sensitivity	92.3%
Specificity	90.5%
Reagents	<ul style="list-style-type: none"> • Membrane test assembly impregnated with colloidal gold conjugated to recombinant <i>Helicobacter pylori</i> antigen (tracer) and rabbit IgG • Recombinant <i>Helicobacter pylori</i> antigen (capture) and anti rabbit antisera at the respective regions
Test duration	Fifteen minutes assay
Storage / Stability	4-30°C, 2 years

Presentation	Pack	Cat.No.
Helicocheck kit contains individually pouched device with sample dispensing droppers	10T	501070010
	25T	501070025