



	<h2>Rapid Immunochromatographic Test for Hepatitis A infection</h2>
Application	For the detection of IgM antibodies to HAV in human serum / plasma
Principle	Self performing, two site sandwich immunoassay, lateral flow device
Sensitivity	100%
Reagents	<p>A. Each INSIGHT HAV-IgM kit contains individual pouches each containing a</p> <ol style="list-style-type: none"> 1. Device: Membrane test assembly impregnated with HAV antigen colloidal gold conjugate, mouse IgG colloidal gold conjugate, anti human IgM antibodies and anti mouse IgG antibodies at the respective regions. 2. Desiccant pouch. 3. Sample dropper. <p>B. Sample running buffer.</p> <p>C. Package insert.</p>
Test duration	15 mins.
Storage / Stability	4-30°C, 18 months

Presentation	Pack	Cat.No.
Insight™ HAV-IgM kit contains individually pouched device	10 Tests	10501010