



Dilute Russell's Viper Venom Time (dRVVT) based Lupus Anticoagulant detection system

Application	For screening and confirmation of Lupus Anticoagulants in plasma	
Sensitivity	Phospholipid dependant Lupus Anticoagulants	
Specificity	Not affected by contact and antihaemophilic factor deficiencies or corresponding inhibitors	
Reagents	• Screening Reagent	Phospholipids Lyophilised with dRVV
	• Confirmation Reagent	Enriched Phospholipids Lyophilised with dRVV
	• Calcium Chloride	Ready to use 0.025M solution, insensitive to heparin upto 0.4 IU/ml
Normal Range	<ul style="list-style-type: none"> • With Screening Reagent 28 - 45 seconds • With Confirmation Reagent, 28 - 40 seconds • Results may also be expressed as Ratio 	
Storage / Stability	2-8°C, 12 months	
In Use Stability	Reconstituted stability 2-8°C, 5 days for screening and confirmation reagent	

Presentation	Pack	Cat.No.
LA Screen	5 Tests	10672005

