



	ABO & Rho (D) Forward grouping confirmation card
Application	For confirmation of ABO and Rho (D) blood group systems
Card Configuration	A-B-D-A-B-D / Two profiles per card
Reagent	Matrix™ ABO & Rho (D) forward grouping confirmation card contains: <ul style="list-style-type: none"> • Six microtubes prefilled with Matrix™ gel in suitable buffer, and with specific monoclonal Anti-A, Anti-B and Anti-D (VI-) antibody in appropriate positions
Standardisation	Follows FDA guidelines
Storage / Stability	4-25°C, 12 months

Presentation	Pack	Cat.No.
Matrix™ ABO & Rho (D) forward confirmation grouping card	24 cards	10255024

