



**matrix**<sup>™</sup>  
 GEL•SYSTEM  
**OCTOPLUS**  
 Complete Grouping Card

## Matrix™ Octoplus Complete Grouping Card

<b>Application</b>	For determination of ABO and Rho (D) blood group systems with autocontrol and reverse grouping for confirmation of ABO grouping
<b>Card Configuration</b>	A-B-AB-D-Ctrl-A1-B-O / One profile per card
<b>Reagent</b>	Matrix™ Octoplus Complete Grouping Card contains: <ul style="list-style-type: none"> <li>• Monoclonal Anti-A, Anti-B, Anti-AB and Anti-D (VI-) prefilled with gel from microtube 1-4</li> <li>• Microtube 5 contains neutral gel which serves as a negative control</li> <li>• Microtube 6,7 and 8 contains neutral gel for reverse grouping</li> </ul>
<b>Standardisation</b>	Follows FDA guidelines
<b>Storage / Stability</b>	4-25°C, 12 months

Presentation	Pack	Cat.No.
Matrix™ Octoplus Complete Grouping Card	24 Cards	102820024