



## Latex Agglutination Slide Test

<b>Application</b>	For the detection and semiquantitation of C-Reactive Protein in serum
<b>Sensitivity</b>	~ 0.6 mg/dl
<b>Specificity</b>	Anti-CRP antibodies of IgG class specific to CRP
<b>Latex reagent</b>	Uniform suspension of polystyrene latex particles coated with Anti-CRP antibodies
<b>Positive Control</b>	Reactive with the reagent
<b>Negative Control</b>	Non reactive with the reagent
<b>Test duration</b>	Two minutes assay
<b>Storage / Stability</b>	2-8°C, 18 months

Presentation	Pack	Cat.No.
<b>RHELAX-CRP kit contains reagents and accessories</b>	5 ml	10420005
	25 Tests	10420025
	35 Tests	104200035
	50 Tests	10420050
	70 Tests	104200070
	100 Tests	10420100