



Sensitive Thromboplastin Reagent

Application	For the Prothrombin Time determination in plasma, monitoring Oral Anticoagulant therapy and assessment of Extrinsic Pathway
Reagent	Ready to use Liquid stable, Calcified Thromboplastin reagent
Source	Rabbit Brain
Normal Range	11 - 15 seconds
ISI value	Optimised at ~ 1.0 (1.0 -1.15)
Storage / Stability	<ul style="list-style-type: none"> • 2-8°C, 18 months • 18 - 25°C, 1 week • 37°C, 2 days

Presentation	Pack	Cat.No.
UNIPLASTIN®	8 ml	10620008
	4 ml	10620004
	5 ml	10620005
	12 x 5 ml	10620125