



Virdict4

Rapid Immunochromatographic Test System Pack for Simultaneous Detection of HIV 1/2 antibodies, HCV antibodies, Syphilis antibodies and HBsAg in human Serum/Plasma

Application

Rapid Immunochromatographic test for simultaneous detection of HBsAg and antibodies to HIV, HCV, and Syphilis in human serum/plasma.

Sensitivity

100% for HIV 1 & 2, HCV, HBsAg and Syphilis

Specificity

100% for HIV 1 & 2, HCV, HBsAg and >95% for Syphilis

Reagents

Ready to use individually pouched rapid devices striped with

- **HBsAg:** Anti-HBsAg-colloidal gold conjugate, rabbit IgG-Colloidal gold. Anti-HBsAg at test "T". Goat Anti-rabbit IgG at 'Control' region
- **HIV:** HIV specific recombinant antigen-colloidal gold conjugate, rabbit IgG-Colloidal gold. HIV specific recombinant antigen at region 'T'. Goat Anti-rabbit IgG at 'Control' region
- **HCV:** multi-epitope HCV recombinant peptide antigen- colloidal gold conjugate, rabbit IgG-Colloidal gold. multi-epitope HCV recombinant antigen at 'T'. Goat Anti-rabbit IgG at 'Control' region
- **Syphilis:** recombinant Treponema pallidum antigen-colloidal gold conjugate, rabbit IgG-Colloidal gold. Recombinant Treponema pallidum antigen at 'T'. Goat Anti-rabbit IgG at 'Control' region

Test duration

30 minutes

Storage / Stability

4-30°C, 24 months

Presentation	Pack	Cat.No.
Virdict 4 kit contains,individually pouched test devices, sample dispensing dropper & diluent buffer	25T	402090025
	50T	402090050

