



The Hand-Held Vein Finding Device

KEY FEATURES

- Reduce multiple needle prick & improves site selection for needle insertion.
- No treatment delay due to difficult venous access
- Reduce patient pain due to needle insertion in muscle/ tingling sensation due to hitting of needle in nerve.
- Increase patient comfort & satisfaction. Avoids hematoma.

TECHNICAL SPECIFICATIONS

Number of LED's	22 Orange, 8 Red
Power ON	Indication by Green LED
Low power	Indication by Red LED
Functional opening	17mm
Viewing diameter	33mm
Rechargeable Battery	Lithium Ion, 3.7V, 1150mAh
Battery charger	110-240 volts AC input, 5V DC output 750mA
Battery Backup	1½ hrs. Continuous usage.
Dimensions	L 125 mm; B 65 mm; H 22.7 mm
Weight	90 gms.
Cat No.	825VS000000

LIST OF COMPONENTS SUPPLIED WITH THE INSTRUMENT

Components	Qty.
Disposable Plastic Sleeves	50 Nos
User Manual	1 No.

ACCESSORIES/CONSUMABLES FOR USE WITH VEINSPY

Accessories/Consumables	Cat No.	Pack Size
Disposable Plastic Sleeves	825DS5000050	50 Nos