



## RAL Stainer

RAL Stainer is an automated bath staining unit, compact and efficient. Thanks to its RFID system it eliminates any protocol errors. You can achieve high quality, standardised and reproducible staining through the combination of calibrated RAL stains and a robust reliable mechanism. RAL Stainer is the ideal flexible solution for laboratories to handle their routine analyses.

RAL Stainer is only compatible with RAL Diagnostics staining kits. The recommended staining protocols have been optimised by RAL Diagnostics in collaboration with reference laboratories.

### KEY FEATURES

#### Intelligent

Pre-programmed, customisable protocols

#### Efficient




Low reagent consumption

#### Safe

Sealed staining vessels

#### Simple

Ready-to-use reagents



## RAL Stainer®

### AVAILABLE KITS



Kit RAL Stainer MCDh™  
REF.: 360200-0000 | 500 slides  
Kit RAL Stainer Gram-Nicolle  
REF.: 360220-0000 | 500 slides  
Kit RAL Stainer Gram-Hücker  
REF.: 360210-0000 | 500 slides  
Kit RAL Stainer 555  
REF.: 360250-0000 | 300 slides

### TECHNICAL SPECIFICATIONS

#### DIMENSIONS

- Height: 58,7 cm
- Diameter: 40 cm

#### WEIGHT

- 28 kg
- 61,7 lbs

#### FUNCTIONS

- 2 loading / unloading / drying zones
- 5 stations, (4 staining stations & 1 rinsing station)
- 2x10L containers for rinsing liquid & waste
- Automatic waste draining in 10L
- Number of tests possible: 500 slides
- Time: 11 m
- Life span: 2 weeks after opening

#### SUPPLY VOLTAGE RANGE

110-230V / 50-60HZ / 150 VA

### Authorized Distributor

#### HORIBA India Private Limited

246, Okhla Industrial Estate Phase III, New Delhi - 110020, India, Tel: +91 (11) 4646 5000 / 4669 5001  
<http://www.horiba.com>, Email: [hematology.hin@horiba.com](mailto:hematology.hin@horiba.com)