

VISION Pro

Automatic ESR Analyzer



KEY FEATURES

Innovative Technology

- ☞ Direct testing with EDTA tubes
- ☞ Quick infrared scanning
- ☞ Good correlation with Westergren method
- ☞ Results in 20 minutes

Advanced Features

- ☞ LCD touch screen
- ☞ Dynamic sedimentation curve available
- ☞ Internal barcode reader and thermal printer
- ☞ Easy recharge through dongle

**EASE OF PROCESSING
REDEFINED**



Specifications

Test Speed	20 min/cycle
Model Types	8/16/32 sample positions
Mode	Cycle & Random
Test tube Requirement	EDTA tube with a diameter of 12/13mm
Sample Volume	1.5ml-3.5ml
Display	7 inches LCD touch screen
Barcode Reader	Internal barcode reader
Printer	Built-in Thermal printer
LIS	Bidirectional connect
ESR Dynamic Curve	Yes
Internal QC Management	Yes
Temperature Correction	Yes
Internal Mixing	Yes
Power Supply	100-240V, 50/60Hz
Input Power Supply	220VA
Dimension	352mmx365mmx325mm
Weight	10.5kg

Marketed by:



CPC Diagnostics Pvt. Ltd.

9, 5th Floor, Gokul Tower 9 & 10, CP Ramaswamy Road,
Alwarpet, Chennai - 600 018.

Phone : +91 44 24993989, Email : info@cpcdiagnostics.in

www.cpcdiagnostics.in

Manufactured by:

YHLO