



Proven Reliability, Greater Flexibility

Clinitek Advantus™ Analyzer

Answers for life.

SIEMENS



How Can You Enhance Productivity While Delivering Uncompromised Quality?

With the Clinitek Advantus Urine Chemistry Analyzer

How do you respond to the demand for high productivity and high quality?

With a reliable, high-throughput urine chemistry analyzer that builds on the proven performance of its predecessor—already the established standard for the industry.

How can you ensure consistent high-quality results with an ever-increasing workload?

With a urine chemistry analyzer that automatically prompts technicians to perform regular QC testing, and locks out users until the quality controls pass. That's quality without compromise.

How can you get your technology to work harder, so you don't have to?

With a urine chemistry analyzer that streamlines your workflow with automatic calibration, immediate start-up, a wide range of testing options, flexible reporting, and an intuitive, easy-to-use interface that puts you in complete control of everything.

The Clinitek Advantus Analyzer offers a broad menu of urine tests that help healthcare professionals deliver quality patient care and aid in the diagnosis of numerous diseases and conditions, such as diabetes, kidney diseases/disorders, urinary tract infections, and liver function.





STAT capabilities



Color touch screen



Reliable Quality with High Throughput for Enhanced Productivity

The Clinitek Advantus Analyzer enables automated reading of the Multistix® family of urinalysis test strips in medium to high throughput settings. With throughput of up to 500 tests per hour, the Clinitek Advantus Analyzer delivers all the reliable functionality of its highly respected predecessor, the industry standard Clinitek® 500 Analyzer, and “raises the bar” to an even higher level of testing flexibility.

Microscopy Consolidation

- Automatically flags and prints reports for samples requiring follow-up microscopy
- Technician enters microscopy results via touch screen or PC keyboard
- Consolidated report can be printed and/or transmitted to LIS

Improved Quality Control

- Automatically prompts technicians to perform QC testing
- QC lock-out capability ensures compliance
- Testing schedule based on your quality assurance program
- Automatic calibration

Flexible Operation

- One-touch switching between strip configurations
- Downloads and prints load list from LIS
- Easily accommodates physician requests and STAT testing
- Network ready
- Up to three different user-defined confirmatory flag reports

Proven Reliability & Support

- The most widely used brand of urine analyzers and strips worldwide
- Builds on the proven technology of the Clinitek 500 Analyzer
- Consistently reliable high-quality results clinicians can depend on
- World-class service and support



One-touch switching between configurations



Microscopic results included in printed and LIS reports

The Clinitek Advantus™ Analyzer Tests for:

Leukocyte
Nitrite
Protein
Blood
Glucose
Ketone
Bilirubin
Urobilinogen
pH
Specific Gravity
Creatinine*
Protein-to-Creatinine Ratio*

Siemens Healthcare Diagnostics, a global leader in clinical diagnostics, provides healthcare professionals in hospital, reference, and physician office laboratories and point-of-care settings with the vital information required to accurately diagnose, treat, and monitor patients. Our innovative portfolio of performance-driven solutions and personalized customer care combine to streamline workflow, enhance operational efficiency, and support improved patient outcomes.

Clinitek, Clinitek Advantus, and all associated marks are trademarks of Siemens Healthcare Diagnostics Inc. All other trademarks and brands are the property of their respective owners. Product availability may vary from country to country and is subject to varying regulatory requirements. Please contact your local representative for availability.

Specifications subject to change. Please refer to the product specification sheet for the most up-to-date information.

Global Siemens Headquarters

Siemens AG
Wittelsbacherplatz 2
80333 Muenchen
Germany

**Global Siemens
Healthcare Headquarters**

Siemens AG
Healthcare Sector
Henkestrasse 127
91052 Erlangen, Germany
Phone: +49 9131 84 - 0
www.siemens.com/healthcare

Global Division

Siemens Healthcare Diagnostics Inc.
511 Benedict Avenue
Tarrytown, NY 10591-5005
USA
www.siemens.com/diagnostics

www.siemens.com/diagnostics